

## Quick Guide #6C: ZOLL M Series Configuration and Operation with Rosetta-DS (Serial Cable Connection)

### I. Configuration instructions for ZOLL M Series for use with Rosetta-DS via Serial Cable

**Note:** The M Series monitor must have the “RS-232” serial/ECG port and it is recommended to have software version “38.50” or newer to ensure compatibility. Contact your ZOLL representative for information on obtaining software updates.

- Required settings:
  - Upload Baud Rate - 38400
  - Enable Data Relay<sup>3</sup> - Yes
  - Auto Transmit After 12 Lead – Yes<sup>1</sup>
- Entering System Configuration<sup>2</sup>

After the monitor has been off for 10 seconds or more;

- Simultaneously press and hold the rightmost and leftmost softkeys.
- Turn the selector switch to **MONITOR** or **ON** (softkeys still pressed).
- Continuing to hold the buttons depressed until the **Passcode** screen appears.

- Enter Passcode

Enter the Passcode (typically 00000000).

- Press the 0 (zero) softkey. This will move the highlight to the next digit.
- Repeatedly press the 0 (zero) softkey until all eight digits have been entered.  
**Note:** If the Passcode was changed from the default, enter the new code instead.
- The Configuration Passcode screen will then automatically advance to the **Select Language** menu.
- Press the **Enter** softkey to select the highlighted language (**ENGLISH**) to move to the main configuration menu.

- Changing Settings

- Press the **Change Config** softkey to display the first configuration option menu.  
**Note:** There are only three user configurable features available within each option menu/screen.
- Keep pressing **Next Item** softkey until you reach each of the features listed in the table below. Once the desired line is highlighted, press the

**Change Value** softkey to change the setting to match the values from the table below.

### Feature Setting

Feature	Setting
ENABLE DATA RELAY <sup>3</sup>	Yes
AUTO TRANSMIT AFTER 12 LEAD	Yes <sup>1</sup>
UPLOAD BAUD RATE	38400 bps

- After configuring the required settings, press the **Exit Config** softkey to display the save/exit menu.


### Notes:


1. The AUTO TRANSMIT feature is not fully functional on software versions older than 38.50.
2. Configuration settings may vary with monitor software version.
3. Data Relay was previously called Interlink by ZOLL. Certain versions of software may still call it Interlink. If Data Relay does not appear in the menus, look for Interlink instead.

## II. Connecting ZOLL M Series to PC with Rosetta-DS

1. Connect the M Series monitor to the computer with Rosetta-DS by using a ZOLL serial cable as shown below.

**Note:** This cable has a 3.5mm plug on one end which plugs into the ECG x1000 RS-232 jack on the M Series and a DB9 F plug on the other end.





### III. Configuration instructions for Rosetta-DS for use with ZOLL M Series via serial cable


**Note:** Computer must have an available RS-232 Serial Port.

1. Launch Rosetta-DS. At the main screen of Rosetta-DS click the [Settings] button. A password box may appear. The default password is “GDrosetta”.



(Continued on next page)






2. Next click the [Setup] button.



3. Check the Enable Comm Port checkbox, select the Comm Port (Usually Comm 1) and set the Baudrate to 38400.




**NOTE: If Comm Port 1 does not work, try Comm Port 2.**

4. Click the  button.
5. To complete the step, click the  button, then click .
6. Click the  button at the top right of the Rosetta-DS window, then click .
7. Restart Rosetta-DS by double clicking the shortcut on the desktop.

(Continued on next page)

#### IV. Acquiring and Sending a 12-Lead from the ZOLL M Series to Rosetta-DS


1. Press the **12 Lead** softkey.



2. Press the **Acquire** softkey to start acquiring the 12-Lead.




3. Once you have acquired a 12-Lead, press the **Dial Phone #** button to transmit the 12-Lead to the Rosetta-DS.



(Continued on next page)

4. The transmission indicator will start moving up and down indicating that the 12-Lead is being transmitted to Rosetta-DS.



**See Quick Guide #9: Sending a 12-Lead from Rosetta-DS to Destinations**

**-END PROCEDURE-**