

CAREpoint Workstation

Application Note: 183V035A


Title: Applying the Power On/Off label for Winsystems PC Modules

Overview: This procedure will show you how to apply the Power On/Off label for Winsystems PC Modules.

- 1) Locate the small power button on the front of the PC Module next to the USB connector. Apply a small amount of rubbing alcohol on a paper towel and then clean the area to the right of the Power Button where the label will be placed



- 2) Place the Power On/Off label to the right of the Power Button. The word "POWER" should line up with the Power LED and the word "ON/OFF" should line up with the Power Button.



-END PROCEDURE-