

B.23 Quick Guide #23: Synchronizing Time To Another Workstation

Enter Pin code...

NOTE: If you see **Click the 'i' button for further details. MSG: PIN01** after entering any # codes listed in this document, contact the CAREpoint administrator for the PIN code. If the PIN is unknown, contact General Devices to obtain a temporary one.

1. On the blue, on-screen keypad, type in #-C-O-N-F-I-G (#2-2-6-3-4-4)



2. Click the after  “Select Program:” and choose Program #3: Network

Select Program :

3. In “Time / Date (Source) Workstation:” enter the name of the CAREpoint Workstation (GDCP-XXXX) that you want time to synchronize to.

Time / Date (Source) Workstation :


Note: This name specifies the ID of the workstation that the Time / Date will be synchronized to.

Note: Get the XXXX number by typing #-A-B-O-U-T (#-2-2-6-8-8) on the unit you want to synchronize to. Then close the About window.

4. Click 

5. On the blue, on-screen keypad, type in #-B-Y-E (#-2-9-3)



6. Click the Start menu, and click  CAREpoint

-END PROCEDURE-