

General Devices Rosetta-DS System Compatibility Information

Note: A signed copy is REQUIRED for acceptance of a purchase order.

Monitor/Defibrillators

The following is a list of compatible monitors / features for use with Rosetta-DS.

- 12-Lead ECG/Summary Data:
 - Philips HeartStart MRx – with Rev 9.00 or later and the “B11 12-Lead Transmission, Rosetta-Lt” option is required for serial cable connection to Rosetta devices. For Bluetooth transmission, the “B10 Event Summary Bluetooth” option is required as well as the “B06 12-Lead Transmission, Bluetooth” option. Contact your Philips representative for upgrading or more information.
 - ZOLL M Series with Bluetooth or “RS-232” serial/ECG port & Interlink/Data Relay for serial cable connection, and Rev 38.5 or later software Note: M-CCT model is NOT compatible.
 - ZOLL E Series with Bluetooth or “RS-232” serial/ECG port & Interlink/Data Relay for serial cable connection, and Rev 1.01 or later software.
 - Physio Control LIFEPAK 12 and LIFEPAK 15 – Requires the additional use of Rosetta-Lt and RCA as an interface to Rosetta-DS (due to Physio restrictions). Contact General Devices for more information.

Note: Monitor features/options/upgrades are NOT quoted or purchased through General Devices.

Data Communications Devices (Radio/Cellular)

- Cellular:
 - Cellular broadband data via Aircard, mobile gateway or WiFi hotspot with internet connectivity to PC upon which the Rosetta-DS software will be installed.
- Digital Trunked Radio:
 - Motorola TRS with Integrated Voice & Data (IV&D) or High Performance Data (HPD), with Radio-IP support and network connectivity to GD-Net (web) server.
 - Harris – Open Sky (compatibility pending approval)
 - Mobile or portable radios must be data capable with Serial or USB connection available

Note: Wireless Connectivity is not provided for by General Devices. Contact your wireless or 2-way provider for more information or services.

Computer Devices

- Laptop/notebook with:
 - 1GHz Pentium-III CPU or faster with at least 512MB RAM
 - 800x600 or higher Display Resolution
 - Microsoft Windows XP or Windows 7 with .NET Framework 3.5 SP1 or higher
 - Default user account needs full access (Read, Write and Delete) to the following folders:
 - Windows XP: C:\Documents and Settings\All Users\Application Data\General Devices\ and all subfolders
 - Windows 7: C:\ProgramData\General Devices\ and all subfolders
 - Network:
 - Internet Connection – TCP: 443 to <https://www.gdserver.net>
 - Local Connection (Between Rosetta-DS and e-Net Data Agent applications on the local PC only, not internet) – TCP: 30001-30010
 - Peripherals: Available Serial RS-232, USB or Bluetooth with WIDCOMM/BroadCOM or Toshiba stack

Note: If using a USB-to-Serial adapter, General Devices recommends using one based on the FTDI chipset. There are known compatibility issues with Windows 7 and USB-to-Serial adapters based on other chipsets.

Receiving Devices

- CAREpoint Workstation with Rosetta-DS option
- CAREpoint-FS (e-mail/fax Forwarding Service)

Contact General Devices' Technical Support regarding compatibility questions at: (201) 313-7075 or support@general-devices.com

I have read and understand the above compatibility information:

Name: _____ Title: _____

Signature: _____ Date: _____

A signed copy is REQUIRED for acceptance of purchase order. Please fax/email a completed & signed copy to General Devices.

1948002I

1000 River Street, Ridgefield, NJ 07657
tel 201.313.7075 fax 201.313.5671
www.general-devices.com
sales@general-devices.com

**GENERAL
DEVICES**

