

If they can hear your voice ...

You can send a 12-Lead!

# Rosetta-Lt<sup>TM</sup>

## 12-Lead ECG & Data Translator

Perfect  
for  
STEMI  
programs!



*Fast and Reliable .....*

*Use radio, landline or cellphone\* to send single & 12-Lead ECGs and other physiologic data - works with Philips, ZOLL and Physio-Control monitors.*

### ***Benefits***

Rosetta-Lt eliminates sole reliance on cellular telephone, giving you greater reliability and redundancy. In addition, you can send reports and talk on the same communications equipment without disruptive delays or multiple calls. Best of all, Rosetta-Lt works where fax and cellular data transmissions may not.

### ***Many Important Features***

As a companion to Philips, ZOLL and Physio-Control monitors, Rosetta "captures" 12-Lead reports & other data for either immediate or off-scene transmission. Rosetta connects to radios or compatible cellphones by cable or, using the built-in acoustic coupler, any landline phone. Compact, lightweight and rugged, Rosetta-Lt easily stows in your monitor's carrying case.

### ***Perfect for STEMI programs***

For STEMI patients, time is everything. With fast transmission and only one button to press, Rosetta puts a fully formatted printout into the physician's hands quickly, reliably and *VERY* easily.

### ***Send Single Lead ECGs***

Rosetta-Lt is also a single lead ECG modulator that is compatible with all standard 1-Volt ECG outputs. Digitally controlled and drift-free, it provides exceptionally high quality printouts.

### ***Choice of Integrated or Stand-Alone Receivers***

In the ED, Rosetta reports may be received on the General Devices' CAREpoint<sup>TM</sup> EMS Workstation or the Rosetta-Rx<sup>TM</sup> stand-alone receiver ... whichever best meets your space and budgetary needs. Each gives you the same crisp, easy-to-read laser printout. CAREpoint's automatic forwarding feature allows you to reliably distribute the 12-Lead report wherever it is needed.

To find out more, visit [www.general-devices.com](http://www.general-devices.com)

Telemedicine Systems • EMS Workstations • Communications Matrixes • Data Translators

