



GD e-Bridge™
GoBOX



Proven, Customizable + Affordable Mobile TOOLS TO EMPOWER EMS + MEDICAL TEAMS

GD IS THE TRUSTED LEADER IN EMS-HOSPITAL COMMUNICATIONS AND MOBILE TELEMEDICINE SYSTEMS FOR OVER 30 YEARS, CONNECTING CARE PROVIDERS NATIONWIDE AND HANDLING THOUSANDS OF EMERGENCY CALLS DAILY.

Gobox is a portable, one-box solution for linking on-scene out-of-hospital care and distant providers with GD e-bridge™ mobile telemedicine. GD e-Bridge™ represents a major advance in responses to natural and man-made disasters, multi-casualty incidents, triage, mobile integrated healthcare and community paramedicine. Medical specialists can be "virtually present" to assist in collaboration with field providers for the right care, in the right place, at the right time. The go-box fully integrates with gd's GD CAREpoint™ for complete solutions.

WHY GD?

- Leader in EMS-Hospital communication since 1990
- Most experienced mobile telemedicine solution provider
- Over 800 Hospital implementations
- Engineering and manufacturing in USA
- Exceptional in-house support with 97% renewal rate
- Award winning company and innovation solutions
- FDA registered with listed devices



WHAT IT INCLUDES

- Rugged enclosure
- Panasonic Toughbook CF-33
- Headset or hands-free audio
- Powered by 12VDC or 115VA
- PTZ Video Camera Controls
- Easy-to-use touch-screen control
- Connections for USB, BT and wireless devices
- Dual-carrier mobile LTE router - FirstNet Ready
- Works with variety of devices and peripherals
- Built-in accessory storage for fast, easy set-up
- GD e-Bridge™ Mobile Telemedicine App with location tracking
- Loud Speaker for quality sound in noisy environments



First Gen GD GoBox



Comprehensive + interactive

